

Anti-NGF/bNGF Antibody (AS2886401-00)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | Anti-NGF/bNGF Antibody (AS2886401-00) is a human monoclonal antibody (mAb) targeting NGF/bNGF. It inhibits intracellular calcium (Ca ²⁺) influx and provides sustained analgesic effects during movement in an osteoarthritis (OA) rat model. This antibody is applicable for research in osteoarthritis. |
| Targets(IC50) | Trk receptor |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481