

SGN-PDL1V Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SGN-PDL1V Antibody is an antibody that targets PD-L1. It can be utilized in the synthesis of ADCs for the compound ADCSGN-PDL1V.
Targets(IC50)	ADC Antibody

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481