

STX-100 (Biotinylated)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	STX-100 (Biotinylated) is a biotin-labeled version of STX-100, which is a humanized antibody expressed in HEK293 cells and targets Integrin $\alpha\text{V}\beta\text{6}$ (ITGAV & ITGB6).
Targets(IC50)	Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481