

Anti-CD371/CLEC12A Antibody (H6e7)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-CD371/CLEC12A Antibody (H6e7) is a human-derived monoclonal antibody (mAb) that specifically targets CD371/CLEC12A. The recommended isotype control for this antibody is Human IgG1 kappa [Isotype Control].
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481