

GT-02216

Chemical Properties


CAS No. : 2648407-76-3

Formula: C₂₂H₂₃N₅O₂

Molecular Weight: 389.45

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	GT-02216 is capable of binding allosterically to GCase, thereby enhancing its function. In primary human fibroblasts, it increases GCase activity and decreases Tau accumulation in mutant GBA1 fibroblasts. Additionally, GT-02216 offers protection to hippocampal primary neurons exposed to Tau oligomer in a rat model.
Targets(IC50)	glycosidase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5677 mL	12.8386 mL	25.6772 mL
5 mM	0.5135 mL	2.5677 mL	5.1354 mL
10 mM	0.2568 mL	1.2839 mL	2.5677 mL
50 mM	0.0514 mL	0.2568 mL	0.5135 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481