

BI-1910

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BI-1910 is a human monoclonal antibody (mAb) targeting TNFRSF1B. It enhances T cell activation and induces CD8 T cell-dependent anti-tumor activity. BI-1910 is applicable for research on solid tumors. Recommended isotype control: HumanIgG1lambda, Isotype Control.
Targets(IC50)	TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481