


Pteridic acid A

Chemical Properties

CAS No. :	367510-96-1
Formula:	C ₂₁ H ₃₂ O ₅
Molecular Weight:	364.48
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pteridic acid A is a plant growth promoter with auxin-like phytohormone activity, capable of enhancing root growth in legumes and other plants at a concentration as low as 1 nM (nanomolar).
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7436 mL	13.7182 mL	27.4363 mL
5 mM	0.5487 mL	2.7436 mL	5.4873 mL
10 mM	0.2744 mL	1.3718 mL	2.7436 mL
50 mM	0.0549 mL	0.2744 mL	0.5487 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481