

Anti-FOLH1/PSMA Antibody (J591)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-FOLH1/PSMA Antibody (J591) (MEDI3726 antibody; ADCT-401 antibody) is a humanized monoclonal antibody (mAb) that specifically targets PSMA. This antibody is utilized in prostate cancer research.
Targets(IC50)	PSMA

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481