

Sym-021

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sym-021 is a human monoclonal antibody (mAb) targeting PDCD1/PD-1/CD279. It blocks the binding of PD-L1 and PD-L2 ligands, which stimulates the secretion of interferon IFN- γ and IL-2 and promotes T cell proliferation. Sym-021 shows antitumor activity in PDX mouse models.
Targets(IC50)	PD-1/PD-L1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481