

NN-8828

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	NN-8828 (NNC-114-0005) is a humanized monoclonal antibody inhibitor that targets interleukin-21 (IL-21). It holds potential for research into immune system disorders.
Targets(IC50)	Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481