

SAR-439459

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	SAR-439459 is a human monoclonal antibody (mAb) targeting TGFB1/TGFβ-1. It blocks TGFβ-mediated pSMAD signaling and inhibits T cell and NK cell activity. Additionally, SAR-439459 exhibits antitumor properties.
Targets(IC50)	TGF-beta/Smad

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481