

TQC2731

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	TQC2731 is a fully human monoclonal antibody that targets TSLP. It can be utilized for research on asthma and other related diseases mediated by Th2 inflammation.
Targets(IC50)	Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481