


## 3-Oxoheptadecanoyl-CoA

## Chemical Properties

CAS No. :	34619-89-1
Formula:	C37H64N7O18P3S
Molecular Weight:	1019.93
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Oxoheptadecanoyl-CoA (3-Oxoheptadecanoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9805 mL	4.9023 mL	9.8046 mL
5 mM	0.1961 mL	0.9805 mL	1.9609 mL
10 mM	0.098 mL	0.4902 mL	0.9805 mL
50 mM	0.0196 mL	0.098 mL	0.1961 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481