

AGEN1307

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	AGEN1307 is a human monoclonal antibody (mAb) targeting TIGIT. It features an active FcγR binding region, enhancing antibody-dependent cellular cytotoxicity (ADCC). AGEN1307 is utilized in cancer research.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481