

TGTB227

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TGTB227 is a human monoclonal antibody (mAb) targeting TIGIT. It works by reducing Foxp3 expression in regulatory T cells (Tregs), thereby inhibiting their immunosuppressive function. TGTB227 is applicable for research in tumor immunology.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481