

AR-201

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	AR-201 is a human IgG1 monoclonal antibody that targets the F glycoprotein F0. It exhibits activity against the respiratory syncytial virus (RSV) and can be utilized in studies related to lower respiratory tract infections (LRTI). The recommended isotype control is HumanIgG1kappa, Isotype Control.
Targets(IC50)	RSV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481