

Anti-HLA-A2 (Hu8F4)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-HLA-A2 (Hu8F4) is a human-derived monoclonal antibody (mAb) that targets HLA-A2. It is applicable in the study of acute myeloid leukemia. The recommended isotype control is Human IgG1 kappa (Isotype Control).
Targets(IC50)	MHC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481