

ADI-65534

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ADI-65534 is a human IgG1 monoclonal antibody targeting Envelopment polyprotein and M polyprotein. It is utilized in hantavirus infection research. The recommended isotype control for this antibody is HumanIgG1kappa [Isotype Control].
Targets(IC50)	Influenza Virus

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481