


DH03

## Chemical Properties

CAS No. : 361162-95-0  
 Formula: C<sub>18</sub>H<sub>14</sub>ClFN<sub>4</sub>O  
 Molecular Weight: 356.78  
 Storage: Keep away from direct sunlight  
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	DH03 is a ligand for target proteins in PROTACs (Ligands for Target Protein for PROTACs). It serves as an active control for the target protein ligand of PROTAC SARS-CoV-2 Degradar-1.
Targets(IC50)	Ligands for Target Protein for PROTAC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8028 mL	14.0142 mL	28.0285 mL
5 mM	0.5606 mL	2.8028 mL	5.6057 mL
10 mM	0.2803 mL	1.4014 mL	2.8028 mL
50 mM	0.0561 mL	0.2803 mL	0.5606 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481