

CDX-1140

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	CDX-1140 is a human IgG2 CD40 agonist antibody. It activates dendritic (DC) and B cells and induces NF- $\kappa$ B stimulation in CD40-expressing reporter cell lines. CDX-1140 functions independently of Fc receptor binding and works synergistically with natural CD40L signaling. It exhibits both direct and immune-mediated antitumor activity, making it useful for lymphoma and solid tumor research.
Targets(IC50)	TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481