


O-Benzyl posaconazole

Chemical Properties

CAS No. :	170985-86-1
Formula:	C44H48F2N8O4
Molecular Weight:	790.90
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	O-Benzyl posaconazole [4-[4-[4-[4-[[[(3r,5r)-5-(2,4-Difluorophenyl)-5-(1,2,4-triazol-1-ylmethyl)oxolan-3-yl]methoxy]phenyl]piperazin-1-yl]phenyl]-2-[(2s,3s)-2-phenylmethoxypentan-3-yl]-1,2,4-triazol-3-one] is an intermediate compound utilized in the synthesis of various active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2644 mL	6.3219 mL	12.6438 mL
5 mM	0.2529 mL	1.2644 mL	2.5288 mL
10 mM	0.1264 mL	0.6322 mL	1.2644 mL
50 mM	0.0253 mL	0.1264 mL	0.2529 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481