

FFP-104

## Chemical Properties

CAS No. : 2378693-56-0

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	FFP-104 is a humanized IgG4 monoclonal antibody inhibitor that targets TNFRSF5/CD40. VX70 is applicable for research into chronic autoimmune liver diseases, such as primary biliary cholangitis (PBC).
Targets(IC50)	TNF

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481