

MEDI-552

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	MEDI-552 is an IgG monoclonal antibody targeting CD20/transmembrane 4 superfamily member A 1 (MS4A1).
Targets(IC50)	Immunology/Inflammation related

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481