


16:0 Folate cap PE disodium

Chemical Properties

CAS No. :	2456348-69-7
Formula:	C62H100N9Na2O14P
Molecular Weight:	1272.46
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	16:0 FolatecapPE is a lipid layer incorporating a phospholipid surfactant. It is applicable for the preparation of lipid nanoparticles used in drug delivery.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7859 mL	3.9294 mL	7.8588 mL
5 mM	0.1572 mL	0.7859 mL	1.5718 mL
10 mM	0.0786 mL	0.3929 mL	0.7859 mL
50 mM	0.0157 mL	0.0786 mL	0.1572 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481