


1-([1,1'-Biphenyl]-2-yl)piperazine hydrochloride

Chemical Properties

CAS No. :	1386928-18-2
Formula:	C ₁₆ H ₁₉ ClN ₂
Molecular Weight:	274.79
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-([1,1'-Biphenyl]-2-yl)piperazine hydrochloride is a pharmaceutical intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6391 mL	18.1957 mL	36.3914 mL
5 mM	0.7278 mL	3.6391 mL	7.2783 mL
10 mM	0.3639 mL	1.8196 mL	3.6391 mL
50 mM	0.0728 mL	0.3639 mL	0.7278 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481