

SHR-1703

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	SHR-1703 is a long-acting humanized anti-IL-5 IgG1 monoclonal antibody. It selectively identifies and binds to IL-5, preventing IL-5 from associating with the $\alpha$ chain of eosinophil surface receptors. This action inhibits the proliferation, activation, and migration of eosinophils, ultimately reducing eosinophil-mediated inflammation and damage. SHR-1703 is applicable for research into eosinophil-related diseases, such as eosinophilic asthma. The recommended isotype control is HumanIgG1kappa, Isotype Control.
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481