

ST2485

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ST2485 is a monoclonal anti-tenascin antibody. It demonstrates additive binding to tenascin C both in vitro and in xenotransplant models. ST2485 binds to human tenascin, sharing an epitope with the BC2 region. It cross-reacts with mouse tenascin and can be utilized in anti-tumor research.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481