

NI-1201

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	NI-1201 is a monoclonal antibody that targets a murine epitope recognized by 25F10 within the human IL-6R sequence. It can inhibit both cis and trans signaling of IL-6. Additionally, NI-1201 specifically targets the IIb site of hIL-6R and prevents the binding of gp130 to IL-6R.
Targets(IC50)	Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481