


Photoswitchable sphingosine

Chemical Properties

CAS No. :	2456348-60-8
Formula:	C ₂₂ H ₂₉ N ₃ O ₂
Molecular Weight:	367.49
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Photoswitchable Sphingosine is a photoswitchable lipid used to study the reversible cellular responses initiated by light stimulation.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7212 mL	13.6058 mL	27.2116 mL
5 mM	0.5442 mL	2.7212 mL	5.4423 mL
10 mM	0.2721 mL	1.3606 mL	2.7212 mL
50 mM	0.0544 mL	0.2721 mL	0.5442 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481