


## (R,S,S)-VH032-Me-glycine

## Chemical Properties

CAS No. :	2407458-67-5
Formula:	C <sub>25</sub> H <sub>35</sub> N <sub>5</sub> O <sub>4</sub> S
Molecular Weight:	501.64
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	(R,S,S)-VH032-Me-glycine is an E3 ubiquitin ligase ligand-linker conjugate, which can be utilized in the synthesis of XF067-68.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9935 mL	9.9673 mL	19.9346 mL
5 mM	0.3987 mL	1.9935 mL	3.9869 mL
10 mM	0.1993 mL	0.9967 mL	1.9935 mL
50 mM	0.0399 mL	0.1993 mL	0.3987 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481