

Anti-IL-6 Antibody (Chugai SK2)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-IL-6 Antibody (Chugai SK2) is a murine monoclonal antibody targeting IL-6. It is applicable in research related to inflammatory diseases and cancer, including rheumatoid arthritis (RA), Crohn's disease (CD), as well as pancreatic, lung, and colon cancer.
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481