

## Pomalidomide-PEG6-C2-COOH

## Chemical Properties

CAS No. : 2225148-49-0


Formula: C<sub>28</sub>H<sub>39</sub>N<sub>3</sub>O<sub>12</sub>

Molecular Weight: 609.62

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Pomalidomide-PEG6-C2-COOH is a synthetic E3 ligase ligand-linker conjugate, consisting of a cereblon ligand based on Pomalidomide and a linker with 6 units of PEG.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6404 mL	8.2018 mL	16.4037 mL
5 mM	0.3281 mL	1.6404 mL	3.2807 mL
10 mM	0.164 mL	0.8202 mL	1.6404 mL
50 mM	0.0328 mL	0.164 mL	0.3281 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481