

trans-Glutaconic acid

Chemical Properties


CAS No. : 628-48-8

Formula: C₅H₆O₄

Molecular Weight: 130.10

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	trans-Glutaconic acid (trans-2-Pentenedioic acid) is the trans isomer of Glutaconic acid, an organic compound featuring two carboxylic acid groups.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.6864 mL	38.432 mL	76.864 mL
5 mM	1.5373 mL	7.6864 mL	15.3728 mL
10 mM	0.7686 mL	3.8432 mL	7.6864 mL
50 mM	0.1537 mL	0.7686 mL	1.5373 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481