

HS-20093 Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	HS-20093 Antibody is an antibody-drug conjugate (ADC) targeting B7-H3, demonstrating antitumor activity. It is applicable for research involving small cell lung cancer, relapsed or refractory osteosarcoma, and advanced solid tumors.
Targets(IC50)	Antibody-Drug Conjugates (ADCs),Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481