

## Anti-IL6 Antibody (MH166)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-IL6 Antibody (MH166) is a mouse-derived monoclonal antibody targeting human IL-6. This antibody can fully neutralize human IL-6 activity in vitro. However, when combined with hIL-6 in vivo, it enhances hIL-6 activity, extends its half-life, and increases the production of antigen-specific antibodies.
Targets(IC50)	Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481