


## Cronexitide lanocianine

## Chemical Properties

CAS No. :	2041574-23-4
Formula:	C73H92N14O15S2
Molecular Weight:	1469.73
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Cronexitide lanocianine is a diagnostic imaging agent.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6804 mL	3.402 mL	6.804 mL
5 mM	0.1361 mL	0.6804 mL	1.3608 mL
10 mM	0.068 mL	0.3402 mL	0.6804 mL
50 mM	0.0136 mL	0.068 mL	0.1361 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481