


## Clomoxir-CoA

## Chemical Properties

CAS No. :	81489-57-8
Formula:	C <sub>35</sub> H <sub>51</sub> ClN <sub>7</sub> O <sub>18</sub> P <sub>3</sub> S
Molecular Weight:	1018.26
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Clomoxir-CoA (Clomoxir coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9821 mL	4.9103 mL	9.8207 mL
5 mM	0.1964 mL	0.9821 mL	1.9641 mL
10 mM	0.0982 mL	0.491 mL	0.9821 mL
50 mM	0.0196 mL	0.0982 mL	0.1964 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481