


## D2-malonyl-CoA

## Chemical Properties

CAS No. :	60051-91-4
Formula:	C <sub>24</sub> H <sub>38</sub> N <sub>7</sub> O <sub>19</sub> P <sub>3</sub> S
Molecular Weight:	855.59
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	D2-malonyl-CoA (D2-malonyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1688 mL	5.8439 mL	11.6878 mL
5 mM	0.2338 mL	1.1688 mL	2.3376 mL
10 mM	0.1169 mL	0.5844 mL	1.1688 mL
50 mM	0.0234 mL	0.1169 mL	0.2338 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481