


Cholyl-CoA

Chemical Properties

CAS No. :	6050-48-2
Formula:	C45H74N7O20P3S
Molecular Weight:	1158.09
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Cholyl-CoA (Cholyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8635 mL	4.3175 mL	8.6349 mL
5 mM	0.1727 mL	0.8635 mL	1.727 mL
10 mM	0.0863 mL	0.4317 mL	0.8635 mL
50 mM	0.0173 mL	0.0863 mL	0.1727 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481