

PF-04840082

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PF-04840082 is a humanized IgG2 monoclonal antibody inhibitor targeting Dkk-1. It binds to the physiological antagonist (Dkk-1) of the Wnt/LRP5 signaling pathway, which enhances bone mass and density by activating osteoblasts. PF-04840082 is applicable for osteoporosis research.
Targets(IC50)	Wnt/beta-catenin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481