

## Pentane-1,3,5-triol

## Chemical Properties

CAS No. :	4328-94-3
Formula:	C <sub>5</sub> H <sub>12</sub> O <sub>3</sub>
Molecular Weight:	120.15
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Pentane-1,3,5-triol is a pharmaceutical intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.3229 mL	41.6146 mL	83.2293 mL
5 mM	1.6646 mL	8.3229 mL	16.6459 mL
10 mM	0.8323 mL	4.1615 mL	8.3229 mL
50 mM	0.1665 mL	0.8323 mL	1.6646 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481