

Arcotatug

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Arcotatug is a humanized monoclonal antibody targeting CLDN18.2. Arcotatug is applicable for the synthesis of ADCs (Antibody-Drug Conjugates). |
| Targets(IC50) | ADC Antibody,Tight Junction Protein |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481