


## 3-Methyl-4-(pyridin-4-yl)aniline

## Chemical Properties

CAS No. :	82160-27-8
Formula:	C <sub>12</sub> H <sub>12</sub> N <sub>2</sub>
Molecular Weight:	184.24
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Methyl-4-(pyridin-4-yl)aniline is a pharmaceutical intermediate employed in the synthesis of various active compounds.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4277 mL	27.1385 mL	54.277 mL
5 mM	1.0855 mL	5.4277 mL	10.8554 mL
10 mM	0.5428 mL	2.7139 mL	5.4277 mL
50 mM	0.1086 mL	0.5428 mL	1.0855 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481