

MDX-1411

Chemical Properties

CAS No. : 2101467-47-2

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	MDX-1411 is a fully human monoclonal antibody that targets CD70. It is capable of inducing antibody-dependent cellular cytotoxicity (ADCC) and can be utilized in the study of B cell malignancies.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481