


3'-NH₂-CTP

Chemical Properties

CAS No. : 85708-23-2
 Formula: C₉H₁₇N₄O₁₃P₃
 Molecular Weight: 482.17
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	3'-NH ₂ -CTP is a nucleotide analog modified with an amino group at the 3' position of CTP.
Targets(IC50)	Nucleoside Antimetabolite/Analog

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.074 mL	10.3698 mL	20.7396 mL
5 mM	0.4148 mL	2.074 mL	4.1479 mL
10 mM	0.2074 mL	1.037 mL	2.074 mL
50 mM	0.0415 mL	0.2074 mL	0.4148 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481