


## 2,2-Dibenzyl-1,3-propanediol

## Chemical Properties

CAS No. :	31952-16-6
Formula:	C17H20O2
Molecular Weight:	256.34
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2,2-Dibenzyl-1,3-propanediol is a pharmaceutical intermediate utilized in the synthesis of various active compounds.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9011 mL	19.5053 mL	39.0107 mL
5 mM	0.7802 mL	3.9011 mL	7.8021 mL
10 mM	0.3901 mL	1.9505 mL	3.9011 mL
50 mM	0.078 mL	0.3901 mL	0.7802 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481