

ASP5094

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ASP5094 is a humanized anti-integrin $\alpha 9$ monoclonal antibody. It exhibits excellent safety and tolerability. ASP5094 can be utilized in the research of diseases such as rheumatoid arthritis.
Targets(IC50)	Integrin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481