

PF-04382923

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | PF-04382923 is a humanized IgG2Δa monoclonal antibody. It binds with high affinity to the free C-terminal region of amyloid β (A β) peptides A β 1-40 and A β 1-42. PF-04382923 is utilized in research concerning geographic atrophy (GA) associated with dry age-related macular degeneration (ARMD). |
| Targets(IC50) | Beta Amyloid |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481