

DS-8895

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DS-8895 is a monoclonal antibody targeting EphA2. It enables non-invasive measurement of EphA2 expression within tumors in vivo. DS-8895 can be utilized for studies involving advanced solid cancers that are EphA2 positive.
Targets(IC50)	Ephrin Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481