


## 5-cis-8-cis-Tetradecadienoyl-CoA

## Chemical Properties

CAS No. :	68134-76-9
Formula:	C <sub>35</sub> H <sub>58</sub> N <sub>7</sub> O <sub>17</sub> P <sub>3</sub> S
Molecular Weight:	973.86
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	5-cis-8-cis-Tetradecadienoyl-CoA (Coenzyme A, S-5,8-tetradecadienoate, (Z,Z)-) is a derivative of Coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0268 mL	5.1342 mL	10.2684 mL
5 mM	0.2054 mL	1.0268 mL	2.0537 mL
10 mM	0.1027 mL	0.5134 mL	1.0268 mL
50 mM	0.0205 mL	0.1027 mL	0.2054 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481